

# iziPositioner

The iziPositioner is an in-vehicle posture support seat. It provides stable and comfortable postural support for wheelchair and buggy users who struggle to sit upright independently while being transported in a vehicle. It can be fitted into a range of different vehicles using an adjustable rear strap system. It is a useful addition to a 24 hour positioning routine for the user.

NOTE: The vehicle positioner is not a safety seat, it only provides postural support. It must be used in conjunction with a vehicle's safety belt system.

*“stable and comfortable postural support”*



Height adjustable headrest



Adjustable chest strap



Attachment straps to fit most vehicles



Contoured seat cushion

*“can be fitted into a range of different vehicles”*



**SHONAQUIP®**

## IZIPOSITIONER SPECS

	BABY	SMALL	MEDIUM	LARGE
Back height	480mm	600mm	700mm	800mm
Back width	280mm	320mm	350mm	400mm
Width between Side Supports	140mm	180mm	220mm	250mm
Seat length	280mm	300mm	350mm	400mm

UHAMBO



SHONAQUIP®

**Uhambo Foundation and Shonaquip are a hybrid social enterprise working together towards an inclusive society for children with mobility and other disabilities.**

**Shonaquip's role** is to ensure that children with mobility disabilities have easy access to a range of appropriate assistive device solutions, clinical and technical support services, which ensures a reduction in secondary health complications, improves physical function and the possibility of children being included in everyday activities with their families, their community and society as a whole. Uhambo Foundation is a shareholder of Shonaquip and receives all profits generated.

[info@shonaquip.co.za](mailto:info@shonaquip.co.za) | Cape Town Tel: +27 21 797 8239 | Gauteng Tel: +27 12 665 1211  
[www.uhambofoundation.org.za](http://www.uhambofoundation.org.za) [www.shonaquip.co.za](http://www.shonaquip.co.za)